

Features

- Speak and listen to the operator directly through the pendant
- Easy installation/User friendly
- Operates anywhere from 300 feet to 600 feet* in or around the home †
- Water-resistant – may be worn in bath or shower**
- Can be worn around the neck, clipped to your belt, or around your wrist.
- Battery provides up to six hours of talk time and up to one year of stand-by time.
- 24/7 monitoring with EMT-certified professionals
- Bi-lingual equipment and operators available

*Line of sight.

† Building structure and conditions may affect range.

**Base station should never be exposed to water.

Your MediPendant™ System Includes:

- Two-way voice speakerphone pendant
- Easy to follow set-up guide
- Belt, neck and wrist attachments
- 1 battery backup
- 1 phone connection kit
- 1 power supply



Questions & Answers

Q: What is the advantage of MediPendant™?

A: MediPendant™ enables you to **speak and listen directly through the pendant**. There is no need to be near a base station.

Q: Is MediPendant™ difficult to install?

A: It is very simple to install. Just plug in the power cord and the phone jack. There are no tools required.

Q: Will my phone bill increase?

A: No, the system dials a pre-programmed toll-free number.

Q: What if there is a loss of power?

A: The base station has backup battery power which allows the system to be used if the power fails.

Q: Who answers the call when the button is pressed?

A: An emergency response call center operator will answer the call. All of our operators are EMT-certified and trained to provide on-call medical assistance.

Q: What happens if I am unable to speak?

A: When the button is pressed, you will hear a voice recording letting you know that your call for help has gone through. Within seconds, an operator will come on the line. If the operator cannot hear you, help will be dispatched immediately.

HOW DO I ORDER?

Contact your local dealer.

AFFIX LABEL HERE

#: _____

Rev:TriDir 01

The **MEDICAL ALERT SYSTEM**
that allows you to
SPEAK and LISTEN DIRECTLY
THROUGH THE PENDANT



MediPendant™
Experience the Difference

AFFIX LABEL HERE

